

ohmic instruments co.

508 August Street, Easton, Maryland 21601

(410) 820-5111 (800) 626-7713 Fax (410) 822-9633

www.ohmicinstruments.com

e-mail: ohmic@ohmicinstruments.com

**BIOMEDICAL & ENVIRONMENTAL SENSORS, INSTRUMENTS & CONTROLS SINCE 1969****ULTRASOUND POWER METER MODEL UPM-DT-1AV & UPM-DT-10AV
TIME PROVEN WATTMETERS FOR DIAGNOSTIC & THERAPEUTIC APPLICATIONS**
 **Approved Electronics**

- FIVE POINT CERTIFICATION (NIST TRACEABLE)
- HIGH RESOLUTION: 2 OR 20mW (1.5 or 15 mW in grams mode)
- MEASUREMENT RANGE: 0-30 WATTS
- FREQUENCY RANGE: .5 TO 10MHZ
- MEASURES TOTAL PULSED OR CONTINUOUS POWER
- AUTOMATIC ZEROING & STABILIZATION
- DIGITAL DISPLAY & RS-232 INTERFACE
- EQUIPPED WITH CARRYING CASE

**UPM-DT SERIES WATTMETERS VERIFIED BY NIST & FDA**

OHMIC Instruments Models UPM-DT-1AV & 10AV are precision wattmeters designed to measure the power output of diagnostic and therapeutic ultrasound transducers. The UPM-DT series have been in use worldwide since 1986 and conform to the testing guidelines recommended by FDA, JCAHO, AAMI and AIUM. Power is measured by the radiation method and the accuracy of the UPM-DT's are traceable to NIST.

Principle of Operation

The most reliable and repeatable means of measuring ultrasound power is the radiation force balance method. The transducer to be tested is centered above a 45° air-backed cone target in water media, which is linked to a precision balance capable of resolving 0.1 or 1.0 mg. When acoustic energy is applied to the cone, the resultant force exerted on the load cell of the balance is directly proportional to the total radiated power. The test tank is lined with sonic absorbent rubber to prevent power reflection. The balance is programmed to convert milligram magnitude forces directly to a readout in watts with resolution to the nearest 0.002 watt for Model UPM-DT-1AV and 0.02 watt for Model UPM-DT-10AV. Verification of the UPM-DT-Series is easily

accomplished by placing a calibrated weight on the arm of the target cone. (A 1-gram weight is equal to 14.65 watts.) OHMIC also maintains NIST traceable mass and ultrasound standards, therefore the units can be returned for periodic recertification.

The Importance of Calibrating Ultrasound Transducers

It is important to measure ultrasound power levels to verify correct patient exposure and reduce potential risks. All diagnostic and therapeutic Doppler and ultrasound equipment should be tested periodically. Typical safe values are defined by the temporal average intensity at the point in the acoustic field where the intensity is at its maximum value, referred to as Spatial Peak Temporal Average (SPTA). The FDA and AIUM provide guidelines for the maximum SPTA values in medical applications. Typical safe SPTA values are 10mW/cm² for echo sounding, 100 mW/cm² for Doppler devices and 3W/cm² for therapeutic use. The watt density (watts/cm²) of a transducer is determined by measuring the total power output and dividing by the active cross sectional area of the transducer. A sample form to document and report ultrasound performance is included in the manual.

WARRANTY: All products manufactured by OHMIC Co. are warranted to be free of defects in material and workmanship for one year after delivery. Any equipment found to be defective within this period will be repaired or replaced free of charge.

SPECIFICATIONS - UPM-DT-1AV & 10AV

	UPM DT-1AV	UPM DT-10AV
Power Range	0- 30 Watts	0- 30 Watts
Resolution	2 mW (1.5 mW in grams mode)	20 mW (15 mW in grams mode)
Display Sensitivity	0.002 Watt	0.02 Watt
Accuracy	± 3% @10 Watts, 25°C	± 3% @10 Watts, 25°C
Zeroing Method	Automatic	Automatic
Stabilization	≤3.0 Second Integration	≤3.0 Second Integration
Maximum Transducer Size	3" Diameter	3" Diameter
Transducer Operating Frequency	.5-10MHz	.5-10MHz
Test Media	Degassed Water	Degassed Water
Computer Interface	RS-232, 300-9600 Baud	RS-232, 300-9600 Baud
Default Baud Rate	2400	2400
Power *	* 120 VAC to 12 VAC 400 mA adapter with 6 Ft. cord	* 120 VAC to 12 VAC 500 mA adapter with 6 Ft. cord
EMI Rating	Conforms to CE & FCC	Conforms to CE & FCC
Electrical Safety	Conforms to UL & CSA	Conforms to UL & CSA
Size	11"H x 15" L x 10"W	11"H x 15" L x 10"W
Weight	14 Lbs Net	14 Lbs Net
Carrying Case	Black Anodized Aluminum	Black Anodized Aluminum

* Power packs and cord sets available to meet any configuration or VAC



UPM-DT Series Ultrasound Power Meters feature a high contrast alphanumeric LCD display. Membrane "touch buttons" enable quick setup, automatic zero, selection of units for the power reading (grams or Watts "custom"), and printer control.



The cone target and transducer positioner. (Tank not shown.)



The rear panel of the UPM-DT Series showing RS-232 and printer ports.

ENGINEERING SUPPORT: OHMIC Instruments Co. designs and manufactures a full line of environmental and bio-medical sensors, instruments and controls. Many of our products are custom designed to meet specific requirements. Our engineers will be pleased to discuss your application.

ohmic instruments co.

508 August Street, Easton, Maryland 21601

(410) 820-5111 (800) 626-7713 Fax (410) 822-9633

www.ohmicinstruments.com

e-mail: ohmic@ohmicinstruments.com

